

---

# TeliCam

Intraoral Video Camera

## User Manual

September 1996



# DMTI


DenMed Technologies, Inc.  
1531 W. Orangewood Avenue  
Orange, CA 92868

714-532-1107

Copyright © 1996 by Dental Medical Technology Systems  
*All Rights Reserved*

**TABLE OF CONTENTS**

<b>The TeliCam Intraoral Camera</b>	<b>3</b>
<b>TELICAM INSTALLATION</b>	<b>4</b>
<b>TELICAM OPERATION</b>	<b>6</b>
Powering On	6
Checking and Replacing Fiber Lamp	6
Front Panel Controls	6
Video and Capture Operation	7
Capturing Images	7
Quad Mode	8
Single/Quad Selection from Handpiece	8
<b>HANDPIECE CONTROL</b>	<b>9</b>
Changing Focal Position	9
Extraoral / Arch	9
Intraoral Detent	9
Extreme Close-Up	10
<b>CONTAMINATION AND CLEANING</b>	<b>11</b>
Disposable Sheath Barrier	11
Cleaning	11
<b>CAPTURED IMAGE JITTER</b>	<b>12</b>
Selecting Frame/Field Mode	12
Footswitch Operation	12
<b>HITACHI MODE OPERATION</b>	<b>13</b>
Hitachi Mode Troubleshooting	13
<b>DMTS Technical Support</b>	<b>13</b>
<b>SPECIFICATIONS</b>	<b>14</b>
<b>ACCESSORIES</b>	<b>15</b>
<b>WARRANTY STATEMENT</b>	<b>16</b>



**Figure 1**

TeliCam IntraOral Camera showing handpiece and chassis.


### **The TeliCam Intraoral Camera**

The TeliCam intraoral video camera is a cost-effective high performance dental imaging system with advanced capabilities. Major features include a Teli camera processor, renowned in the intraoral marketplace for excellent color rendition and contrast, a single-lens handpiece not requiring any lens change to view from close-up to extraoral, and integrated video capture memory.

## TELICAM INSTALLATION

Unpack and inspect the TeliCam and any optional equipment received for any signs of shipping damage. Save the packing materials in the event a return must be made. Locate the supplied video and power cables, and become familiar with the connectors on the equipment.

Prepare the location for the TeliCam equipment, and set the units in place. Note that in custom installations adequate clearance must be provided for proper air flow around the TeliCam in order to prevent possible overheating of the fiber light source, as shown in Figure 2. Air is drawn in through the back of the TeliCam chassis, and heated air exhausts out the bottom right side of the unit.



**Figure 2**


Installation inside TeliCam Cart showing proper clearance for air flow.

A thermal switch will shut off the fiber lamp in the event there is an overheating problem. This will be the case if the air flow is blocked or limited.

The TeliCam is a video output device much like a typical camcorder or VCR, and must be hooked to an appropriate video display for images to be seen. The TeliCam has video output connectors for standard composite Video as well as S-Video (Y/C). Select the appropriate output for your installation and hook up the corresponding video cables from the back of the TeliCam to the video display.

Please refer to Figure 3 for a typical cart installation. Be sure to plug the TeliCam into the appropriate building electrical outlet or power strip.

*Note: The TeliCam is intended only to be used with third-wire grounded electrical power. Do not defeat or alter the third-wire safety ground under any circumstances.*



030-1013 PA

**Figure 3**

Typical Cart installation showing TeliCam-Printer-Monitor video cabling.

**TELICAM OPERATION**


**Powering On**

Turn on the TeliCam by pressing “Power” on the left side of the front panel of the unit. A green LED will light up on the power button, and other LEDs will light on the front panel. At this point, the video monitor should show a dark gray screen with a faint black cross in through the center. This is a display of the video memory, which is empty at power up.

Pressing the “Light” button turns on the fiber optic light, and light should be seen emanating from the tip of the handpiece. A corresponding LED will light on the “Light” switch.

**Checking and Replacing Fiber Lamp**

If no fiber light is seen when the “Light” switch is depressed, then it is possible the fiber lamp may have dislodged during shipping. Unplug the TeliCam, and turn it over to reach the lamp door. Remove the lamp door as shown in Figure 4, and check to see that the lamp is fully seated into the socket by gently pressing down on it until resistance is felt. Replace the lamp door and reinstall the system.



**Figure 4**


Bottom of the unit showing the fiber lamp properly installed.

In the event the lamp must be replaced, remove the lamp by actuating the metal arm and sliding the lamp up and out of the socket. Be sure to let the Lamp cool down before touching if it has been in use for a period of time. Refer to the TeliCam Accessories for reorder of the specific lamp type for your model.

**Front Panel Controls**

Many of the features of the TeliCam are controlled via the front panel of the unit, shown in Figure 5. The front panel is made of material which is resistant to most common clinical solutions, and can be easily cleaned by wiping with alcohol or standard cleaning solutions. When activating the ‘Power’ and ‘Light’ switches, a mechanical toggle can be

felt in addition to the visible LEDs. The remaining switches are momentary contact, and actuation is indicated by the corresponding LED.



**Figure 5**

Front Panel of the TeliCam.

The TeliCam is capable of displaying either 'Live' or 'Memory' images on screen. Additionally, the TeliCam can display either a 'Single' image or four 'Quad' images. These modes are selected using the "Live/Memory" and "Single/Quad" toggles located on the right side of the front panel.

### **Video and Capture Operation**

The unit "wakes up" in 'Quad mode' displaying 'Memory'. By pressing one of the capture buttons on the camera handpiece, or the capture pedal, a live video image will be seen on the video monitor. Releasing the button or pedal will freeze the image in one quadrant of the screen. Pressing and releasing the button or pedal will then save images into successive quadrants.

Pressing the "Live/Memory" button on the front panel will toggle the TeliCam between display of the "Memory" image and a "Live" image. By switching to 'Live mode', a live full-screen video image is displayed. Now, after a capture button or pedal is pressed, upon release, the image will freeze for a second or two, and then switch back live. The image has been captured into memory, but live video is still displayed on screen. Selecting 'Memory' from the front panel will then display the captured image on screen.

### **Capturing Images**

There are two "capture" buttons on the TeliCam handpiece. When you are ready to start capturing images, hold hand steady, then press down and lift up on either of the "capture" buttons. This will freeze the picture. The captured picture will quickly return to live, and can be retrieved by pressing "Memory" on the front panel. If this picture is undesirable, simply press down and lift up on the "capture" button again and a new captured image will appear. Pressing "Live" will return the TeliCam to live mode. Images may also be


captured by keeping the TeliCam in "Memory" mode. Keep either handpiece button depressed until ready, then release to capture.

### **Quad Mode**

When 'Quad' mode is selected (default on power-up), successive images will be captured into each of the four quadrants of the screen. A dimly-lit LED will indicate the next quadrant to be captured. When a quadrant is captured, the TeliCam automatically cycles to the next in sequence.

In the event you capture an undesirable picture in one of the quadrants, the four buttons located around the yellow LED on the front panel can be used to select the quadrant you wish to recapture. Press the appropriate "Quad" button. The corresponding LED will light brightly, designating the quadrant which will be recaptured. Retake the new picture until you get it right. Then press the same quadrant button again, and it will toggle back dim indicating the next quadrant to be captured.

An example of the power of the TeliCam Quad mode is shown in Figure 6.



**Figure 6**

An example of four different views displayed simultaneously in Quad mode.

### **Single/Quad Selection from Handpiece**

To change from 'Single' to 'Quad' mode using the handpiece, simply hold one of the "capture" buttons down and press the other at the same time. Repeating this sequence will return TeliCam to 'Single'.



## HANDPIECE CONTROL

### **Changing Focal Position**

There is a movable gold ring on the handpiece which is used to select focal position. Simply pull or push the ring forward or back to set the desired focal position. This can be done with one hand by gripping the hand piece and moving the ring with the thumb or forefinger. The ring will move about 1.5mm. This will allow focusing on objects from as close as 2mm to full face distance, and beyond, in one movement.

### **Extraoral / Arch**

With the ring moved all the way forward (towards the tip of the camera), objects from infinity in to about a full arch or smile will be in sharp focus, as shown in Figure 7.



**Figure 7**

Extraoral / Arch setting (ring forward).

### **Intraoral Detent**

Moving the ring back towards the rear of the handpiece will bring the focus in to closer ranges suitable for intraoral imaging. As the ring is moved back, a slight detent should be felt. This detent is set to a nominal intraoral focal position where 1-2 teeth will be in sharp focus, as shown in Figure 8. Most doctors will find that this position offers a comfortable working position for normal oral examination.




**Figure 8**

Intraoral Setting (ring at detent).

### **Extreme Close-Up**

As the focal ring is moved further back, the range of focus will become closer and closer to the target, and the corresponding image magnification will become higher, approaching x120, as viewed on a 13" monitor. The same molar as in Figure 8 is shown Figure 9, at high magnification.



**Figure 9**

Extreme close-up (ring moved to rear).

**CONTAMINATION AND CLEANING**

**Disposable Sheath Barrier**

The TeliCam is designed to be used in conjunction with a disposable sheath to act as a contamination barrier. Proper installation of the sheath is required to provide the clearest possible image.

When inserting the tip of the camera into the sheath, be sure to keep the camera straight and slide directly into the sheath without twisting. Be sure the camera is aimed down at the paper backing of the sheath packaging in order that the clear window of the sheath line up with the TeliCam lens window. Refer to Figure 10.



**Figure 10**

Properly Installed Disposable Sheath.

**Cleaning**

Care should be taken when cleaning the TeliCam. Do not soak the camera in any solutions, or use any phenol-based disinfectants on the handpiece cable. It is suggested that the handpiece be wiped down with cold sterile or alcohol using a cotton cloth or swab.

If the lens window becomes contaminated, carefully clean with a swab and alcohol. Avoid excessive amounts of any solution to remain on the window area.

## **CAPTURED IMAGE JITTER**

The TeliCam captures single images using both interlaced video fields (“Frame” mode) to ensure maximum capture resolution, by default. If there is movement of the subject during the time these two fields of video information are captured, then the captured image will appear to ‘jitter’ on screen.

To eliminate the possibility of capture jitter, it is possible to set the TeliCam to use only a single field of video (“Field” mode) for capture. Using “Field” mode will eliminate most jitter when capturing a single full-screen picture, but may result in lower resolution than “Frame” mode under certain conditions.

### **Selecting Frame/Field Mode**

Selecting ‘Field’ mode to stop image jitter is done by simply pressing one of the four Quad buttons on the front panel. All four LEDs around the Quad will blink once, indicating one ‘Field’ of video information. Pressing again toggles back to ‘Frame’ mode, the LEDs will blink twice indicating both interlaced video fields are being used.

On some older units, Frame/Field can only be selected after powering the TeliCam on in a special manner:

- With the power off, press and hold down any one of the four Quad buttons on the TeliCam front panel.
- While still pressing a Quad button, power on the TeliCam.

You will then notice all four Quad LEDs blinking until the Quad button is released. The camera is now in “Frame” mode, but pressing any of the Quad buttons while viewing a single image on screen will toggle the image to “Field”, and remove jitter. When toggling from Frame to Field, the Quad LED’s will blink twice to indicate Frame mode, and once for Field mode.

### **Footswitch Operation**

Using the optional footswitch will also help to reduce jitter in captured images by allowing a steadier hand hold on the handpiece. The optional footswitch uses the right pedal to capture images, and the left pedal to toggle from Single to Quad mode. Simply plug the footswitch into the footswitch connector located on the rear of the TeliCam, shown in Figure 3.

### **HITACHI MODE OPERATION**

Provision has been made to tightly integrate the control of the optional Hitachi VY-170A video printer into the TeliCam. By setting the TeliCam to "Hitachi" mode, it is possible to capture in printer memory and initiate printing from the handpiece buttons or foot switch. The Hitachi VY-170A must be linked to the TeliCam via an 8-pin Remote cable, as shown in Figure 3. Simply connect the cable between the TeliCam Remote connector located on the rear of the unit and the Remote connector located at the rear of the printer.

Power up the printer by pressing the "Operate" button on front panel (the green "ready" light will flash at first, then stay on). Press TeliCam 'Live'/'Memory' button for about two seconds, until all four 'Quad' lights turn on. The system is now in "Hitachi" mode.

Press the "Monitor" button on the Hitachi (green light on). Now, the TeliCam handpiece buttons or the footswitch will control the Hitachi printer directly, and bypass TeliCam memory. Capture images by pressing and releasing either hand piece button, or the capture pedal, just as in normal TeliCam operation. Then to print, press both hand piece buttons at once, or the right foot pedal.

To select how many images you would like to appear on the monitor, select the 'Picture' button located behind the Hitachi control panel door. This allows you to split the monitor screen into multiple shots, such as 4 (Quad), or 16.

To leave "Hitachi" mode and return to "Teli Mode", press 'Live'/'Memory' again for about two seconds, until the four 'Quad' lights go out. The TeliCam will now function normally again, and the Hitachi printer will not be affected.

#### **Hitachi Mode Troubleshooting**

If you have followed the above steps and cannot get a picture on the screen:

- Make sure the TeliCam is in "Live" mode (front panel).
- Press the "Input" button on Hitachi, repeat until INPUT S appears.
- Select VIDEO 3 by pressing "Video Select" inside door on front panel of the monitor.
- Double check cable hook-up.

### **DMTS Technical Support**

If you are experience difficulty in this or any other operational issues, please call DMTS Technical Support at (714)-532-1107.

**SPECIFICATIONS**

General

- Intraoral Video Camera NTSC or PAL format
- Integrated Video Freeze Memory with Single/Quad view
- Integrated High Intensity Fiber Optic Light Source
- Handpiece or Footswitch Freeze Memory operation
- Spill resistant membrane switch Front Panel
- Direct Control of Hitachi VY-170A Printer via Remote Connector

Optical

- Focus Range from ~2mm to Infinity
- Intraoral High-Magnification (>120x on 13" monitor)
- Extraoral Full Face capability
- Optical Distortion <±0.5%
- <10mm Dia. Intraoral Tip
- >90mm Intraoral Reach

Video

- |           |                |  |
|-----------|----------------|--|
| • Video   | Composite      | 1 Vpp, 75Ω (unbalanced), SYNC negative |
| • S-Video | Y(Luminance)   | 1 Vpp, 75Ω                             |
|           | C(Chrominance) | 0.28 Vpp, 75Ω                          |

Electrical Power Requirements

NTSC

- 85-136 VAC, 50-60 Hz, 200W Max

PAL

- 180-264 VAC (PAL), 50Hz, 200W Max

Size and Weight

- Chassis Size: 11"W x 9.75"D x 3.30"H
- Weight: 6.5 LB

Environmental (Operating)

- Temperature 5° C to 35° C (40° F to 95° F)
- Humidity 20% to 80% RH

**ACCESSORIES**

The following TeliCam consumables and typical accessories can be ordered directly from DMTS by calling the toll free number 714-532-1107. Be sure to check with your sales representative for other items not listed.

**Consumables**

Disposable Sheaths	50 pk	
Disposable Sheaths	500 pk	
Hitachi Print Film	50 pk	
Toshiba Print Film	50 pk	
Halogen Lamp	1 pk	GE Type ESD, 120V, 150W

**Accessories**

Footswitch, 2 Pedal

Hitachi Printer Remote

Cord, AC Power, 6', US 3-Prong, F-M

Cord, Extension, 9'

Power Strip, 5 Outlet, 15'

Cable, Remote, 6', 8-Pin MiniDIN

Cable, Remote, 20', 8-Pin MiniDIN

Cable, Remote, 40', 8-Pin MiniDIN

Cable, S-Video, 6'

Cable, S-Video, 20'

Cable, Video, 6', RCA-RCA

Cable, Video, 6', BNC-RCA

Wall Mount, 20" Monitor

Wall Mount, 13" Monitor

Pole Mount, PM1 (Universal)

Pole Mount, PM2

Pole Mount, Spacer

**DENTAL MEDICAL TECHNOLOGY SYSTEMS INC.**

**1531 West Oranewood Avenue  
Orange, CA 92868  
714-532-1107**

**WARRANTY STATEMENT**

DMTS warrants the TeliCam product as manufactured by DMTS against defects in craftsmanship and materials for one year from the date the product was delivered to the original purchaser. Computer related hardware is warranted for 90 days. DMTS, at its discretion, will repair or replace any parts of the products (other than users manuals, software and computer related hardware) that are defective in craftsmanship or materials if the product is proven faulty within the one year warranty, provided DMTS is given notice in a timely manner. DMTS will be responsible for shipping costs of any repair/replacement parts for the first 30 days after receipt of a system. The customer will be responsible for shipping costs for any replacement or repair of defective parts thereafter.

If products are to be returned to DMTS, the original customer must obtain an RMA (Return Merchandise Authorization) number from DMTS via telephone. If a return is made after the first 30 days of delivery, the customer must include the RMA number on the package and the customer must prepay all freight charges to return any product to a DMTS designated location. DMTS will ship replacements freight collect on an exchange basis to the customer. Products will become property of DMTS that are returned to DMTS under this warranty. Only the warranty given by the manufacturer will apply in respect to any system or part thereof not supplied by DMTS, if any.

This limited warranty excludes video light bulbs and does not cover damage which is due to: 1) any defects in any product or part not supplied by DMTS as part of the system; 2) fine tuning, relocation, modification or repair by other than DMTS authorized personnel; 3) improper maintenance or installation, carelessness, or negligence or any cause other than ordinary applications, industrial, commercial or personal; 4) any statements made about DMTS products by regional offices, dealers, distributors, agents or salespeople unless confirmed in writing by a DMTS Officer, or; 5) damage which occurs in shipment to or from Customer.

DMTSs' liability shall not exceed price of the individual piece whose loss or defect is the basis of the claim. DMTS is not liable for loss of use of facility, loss of profits, loss of use of equipment or other incidental consequential damages. DMTS is free of obligation beyond its standard warranty with regards to changes made in system component hardware by the manufacturer of that hardware. Without prior notice DMTS reserves the right to make changes in its system hardware at any time. This limited warranty shall constitute the entire limited warranty and is granted by DMTS instead of any or all or all other written, oral, implied or statutory warranties.